

PROPEX[®] II

NEW GENERATION APEX LOCATOR WITH FULL COLOURED DISPLAY



ENDO4YOU

SHARPEN YOUR ROOT CANAL VISIBILITY

MAILLEFER

DENTSPLY

PROPEX® II



PRODUCT OVERVIEW

- Large coloured display to facilitate file tracking
- Progressive sound and visual checks for dual control
- Powered by a rechargeable battery for your convenience
- Latest multi-frequency technology enables location of the apex in most types of conditions


ORDERING INFORMATION

ProPex II	A102800000022
Hooks	A102900000400
Lip clips	A102900000300
Forks	A102900000600
Measuring wire	A102900000500
Battery charger	A102900000222
Battery	A102900000100

DENTSPLY (Australia) Pty Ltd ABN 15 004 290 322
11–21 Gilby Road, Mount Waverley VIC 3149
Tel: 1300 55 29 29 • Fax: 1300 55 31 31
clientservices@dentsply.com • www.dentsply.com.au


UPPER CANAL ZONE

The cursor on the tooth icon indicates the progression of the file in the canal




APICAL ZONE


Extended apical zone from 0.9 to 0.0



When apex is reached, a solid tone is emitted



When the cursor on the tooth icon displays the word 'over', this indicates that the file has passed the apex



DENTSPLY (New Zealand) Ltd
58 Richard Pearse Drive, Airport Oaks, Auckland
Tel: 0800 809 074 • Fax: 0800 809 075
clientservicesnz@dentsply.com • www.dentsply.co.nz