

HU-FRIEDY IZSTĀDES AKCIJAS

PIEDĀVĀJUMS SPĒKĀ LĪDZ 2019. GADA APRĪĻA BEIGĀM!





How the best perform



NEW BLACK LINE SURGICAL COLLECTION

DESIGNED FOR PERFORMANCE





The New Black Line
Surgical Collection completes
the Hu-Friedy solutions for
implant and periodontal
surgeries, with a wider choice
to perform atraumatic extractions,
tunneling techniques and
microsurgical procedures.

THREE MAIN FEATURES BRINGING CRUCIAL CLINICAL SIGNIFICANCE

- Performance engineered coating for enriched contrast, visual acuity and lubricity, reduced light reflection and optimal edge retention
- Matte finished handle for reduced peripheral light reflection and increased clinician comfort
- Unique, smooth, large diameter handle for reduced muscle strain and increased tactile sense

ATRAUMATIC EXTRACTION



EPTSSPX Luxating Hybrid Elevator, Straight Spade

Sugg. Retail Price € 139,00



PT2X
Periotome, Anterior,
Double-Ended

Sugg. Retail Price € 65,00



IMPLANT SURGERY



CRMXMinnesota Retractor

IMP6578SX Kramer-Nevins Sinus Lift, Obtuse

Sugg. Retail Price € 53,50 Sugg. Retail Price € 108,00





NEW BLACK LINE SURGICAL COLLECTION

DESIGNED FOR PERFORMANCE.

TUNNELLING TECHNIQUE



TKN1X Tunneling Knife 1 TKN2X

Tunneling Knife 1

Sugg. Retail Price € 107,00



FQUINCEY1X
Giles de Quincey 1
FQUINCEY2X
Giles de Quincey 2

Sugg. Retail Price € 106,00

SPVXCastro Microsurgical,
Super-Cut

MICRO SURGERY _

Sugg. Retail Price € 523,50



KPAXAllen End-Cutting
Intrasulcular Knife

Sugg. Retail Price € 80,00



PERIODONTAL SURGERY



PPAELAXAllen Periosteal
Elevator, Anterior

Sugg. Retail Price € 110,00



PFIWDS1MKX

Urban Mini Me Periosteal

> Sugg. Retail Price € 48,00



Line instruments.
In addition to
their perfect
quality, they are
supercool.

- DR. ISTVAN URBAN











HU-FRIEDY KEY OPINION LEADER

DR. DOMENICO MASSIRONI MASSIRONI KIT FOR PROSTHETICS BY HU-FRIEDY

Dr. Domenico Massironi, born in Melegnano (MI), graduated with honors in Medicine and Surgery at the University of Pavia. where he subsequently obtained the Specialization with honors in Odontostomatology. He specialized in prosthetic dentistry and has been using the stereomicroscope since 1988, Active member of the EAED (European Academy Esthetics Dentistry) and of AARD (American Academy Restorative Dentistry), he is the author of several publications on prosthetics and precision restoration. Together with his friend and colleague Dr. Carlo Ghezzi, he founded the MEG* training

Visit the website: www.meg-educational.com*

center (Master Educational Group).



The increasing number of aesthetic and prosthesis's patients drives the clinicians to have a wide knowledge and familiarity with techniques and tools that modern dentistry already has, in order to achieve optimal results for their patients and creating a state of wellness that is the aim of aesthetics. The realization of this KIT for prosthetic procedures and of the two new instruments, thanks also to the support of Hu-Friedy, is the demonstration that the best aesthetic results can be obtained for our patients.

-DR. DOMENICO MASSIRONI



The new **MASSIRONI KIT** by Hu-Friedy is set of instrumentation used during prosthetic procedures. Conceived by **Dr. Domenico Massironi**, pioneer of the use of the operating microscope in prosthesis, the kit is designed to be used alongside magnification tools, allowing to finalize prosthetic work with maximum precision and accuracy.

The kit consists of 10 Hu-Friedy instruments, including a tissue plier, a high definition mirror as well as a Klemmer plier, which provide a dynamic feature in modern prosthetic dentistry.

The Massironi Kit for prosthetics includes **2 new Enamel Rounded Chisel** instruments, made with XTS technology, and designed by Dr. Massironi himself in collaboration with Hu-Friedy. These chisels are most often used in the finishing stages of the procedure and the **off-angle design** allows clinicians to access the restoration from any angle.

THE ENAMEL ROUNDED CHISELS BY MASSIRONI

TNMASS1: 1.2 mm chisel - the slender width of this chisel makes it ideal for use during highly magnified prosthodontic cases and minimally invasive preparation techniques.



TNMASS2: 1.5 mm chisel - the slightly larger width of this chisel makes it excellent for everyday use and allows the traditional prosthesis to receive an optimal finishing touch.



MASSIRONI KIT FOR PROSTHETICS BY HU-FRIFDY HOW TO DO?



TNMASS1 & TNMASS2





ENAMEL ROUNDED CHISEL MASSIRONI 1.2 mm & 1.5 mm

The two chisels are most often used in the finishing stages of the procedure and the **off-angle design** allows clinicians to access the restoration from any angle.

PCORTLINI





CORTELLINI, SINGLE END LASER PROBE, .5-15

Ideal for evaluation of minimum thicknesses in the prosthetic reduction of crowns and laminate veneers, and also for the use of silicone index.

GCPYD1





GINGIVAL CORD PACKER YARDLEY #1 HANDLE #41

Spatula for inserting rear deflection wires, also of large diameter, with different angle for optimal work over 360 degrees.

GCPBN1





GINGIVAL CORD PACKER #BN1 HANDLE #41

Spatula for inserting front deflection wires, small-medium diameter, with different angle for optimal work over 360 degrees.

CVIPCKOT4



CARVER KOTSCHY INTERPROXIMAL DOUBLE END #4

Multifunctional ultra-thin spatula, suitable for insertion in the deflection phase of thin silk-like threads and for managing the resin in the margin during the relining of the provisional.



MASSIRONI KIT FOR PROSTHETICS



MIR5HDE

HIGH DEFINITION MIRROR, EU SIMPLE STEM, #5

High definition mirror ideal when using magnifying instruments provided by LED or Xenon lights areas.


МНЕ6СН

EU MIRROR HANDLE SMOOTH SATIN STEEL #6

Used for indirect vision, retraction of lips, cheeks and tongue as well as for illumination of dark areas.

MASSIRONI KIT FOR PROSTHETICS:

PARTCODE	DESCRIPTION	Sugg. Retail Price
MASSKITDIN	MASSIRONI KIT FOR PROSTHETICS	
ARTIKEL-CODE	DAS MASSIRONI-KIT BESTEHT AUS:	
EXD16	Double End Explorer #1, Satin Steel Handle #6	€ 25,00
PCORTLINI	Cortellini, Single End Laser probe, .5-15	€ 27,00
CVIPCKOT4	Carver Kotschy Interproximal Double End #4	€ 35,00
TNMASS1	Enamel Rounded Chisel Massironi, 1.2 mm	
TNMASS2	Enamel Rounded Chisel Massironi, 1.5 mm	
GCPYD1	Gingival Cord Packer Yardley #1 Handle #41	€ 40,00
GCPBN1	Gingival Cord Packer #BN1 Handle #41	€ 40,00
TP32	Semkin-Taylor Tissue Pliers, Curved #32	€ 27,00
Н3	Hemostat Halsted-Mosquito, 12 cm Curved #3	€ 30,00
мне6СН	EU Mirror Handle Smooth Satin Steel #6	€ 21,00
MIR5HDE	High Definition Mirror, EU Simple Stem, #5	€ 12,00
IME12D17C2	IMS Infinity Cassette 1/2 DIN, Black	€ 127,00



XTS TECHNOLOGY

A dark Aluminum Titanium Nitride (AITIN) coating provides a smooth, hard surface which holds a sharp cutting edge. The dark coating, XTS®, also enhances visibility by providing high contrast and also reduces light reflection from operatory and microscope lights.

HU-FRIEDY ADHESTHETICS THECA KIT

BY FFDFRICO FFRRARIS

FEDERICO FERRARIS ADHESTHETICS THECA KIT BY HU-FRIEDY

The Adhesthetics Theca kit is comprised of four unique instruments designed by Dr. Federico Ferraris in collaboration with Hu-Friedy. This kit covers a variety of clinical procedures for both direct and indirect restorative dentistry. Helping clinicians achieve necessary results through clinical excellence. Each instrument has its own distinct set of colour cones, making instrument identification easy. Adhesthetics is the fusion of Adhesion and Esthetics and represents the process by which procedures are codified during the different phases of the treatment. This philosophy gives successful protocols to those clinicians who want to achieve excellence in Restorative Dentistry.



The Kit includes:	CODE	Sugg. Retail Price
Ferraris Main Shaper	TNFF1/2	€ 60,50
Ferraris Fine Shaper/ Cord-packer	TNFF3/4	€ 55,00
Ferraris 3D Shaper	TNFF5/6	€ 55,00
Ferraris Direct Caliper	TNFF7/8	€ 60,50

THE KIT INCLUDES:

- IMS Infinity Cassette 1/4 DIN, IME14DIN58, blue

- 4 Ferraris instruments

TNFFDIN

CODE

Sugg. Retail Price € 331.64









Dr. Federico **Ferraris**

President Elect of SCAD for the biennium 2021-2022

President Elect of AIC for the biennium 2022-2023



adhesthetics:

Incorporating these techniques into everyday procedures allows for predicable results in direct and indirect restorative dentistry. Constantly challenging yourself and continuously making personal improvements is the personal motto of Federico Ferraris.

Do Better Than Yesterday





CHU'S AESTHETIC GAUGES





THE PROPORTION GAUGE

Precise color-coded measurements

- · Provides quick, accurate diagnosis of tooth proportion
- Provides accurate results and reduces chairside adjustment time
- Easy to read reduces visual fatigue

Common reference guide between clinicians and labs

 Results in effective communication to reduce the incidence of errors, and repeated adjustments


Compatible with IMS® cassettes and can be easily sterilized along with other instruments

· Reduces incidence of cross-infection



The T-bar tip has a vertical and horizontal arm. This measures length and width at the same time.







HU-FRIEDY.COM/ChusGauges













Precise color-coded measurements

THE CROWN LENGTHENING GAUGE

- Provides guick, accurate measurements and better results
- Easy to read reduces visual fatigue

Compatible with IMS cassettes and can be easily sterilized along with other instruments

· Reduces incidence of cross-infection



The BLPG tip has predefined measurements to help achieve the proper mid-facial clinical crown and biologic crown length during a crown lengthening procedure.

BLPG TIP



The Papilla tip has predefined measurements to help establish the correct aesthetic position of the interdental papilla from the incisal edge before the flap is closed and sutured.

PAPILLA TIP







PART CODES

PROCLHDLS

PROGS

CLGS SOUNDGS SCHUSET



THE SOUNDING GAUGE

- · Bone sounding made simple and quick
- Sounding tip curvature and sharpness allows easy manipulation and access into deeper areas to analyze the level of the bone crest



The Sounding tip helps determine the sulcus depth. mid-facial osseous crest location and inter-proximal osseous crest location.







DESCRIPTIONS	Sugg. Retail Price
Proportion Gauge (1 Satin Steel Handle, 2 T-Bar Tips, 2 Inline Tips)*	€ 66,00
Crown Lengthening Gauge (1 Satin Steel Handle, 2 BLPG Tips, 2 Papilla Tips)*	€ 66,00
Sounding Gauge (Satin Steel Handle)*	€ 32,00
1 Proportion Gauge, 1 Crown Lengthening Gauge, 1 Sounding Gauge, 15-Instrument IMS Cassette*	€ 210,00
Proportion and Crown Lengthening Gauge Handle (Satin Steel)*	€ 66.00

*Also available in Resin Handle: PROG | CLG | SOUNDG | CHUSET | PROCLHDL

TBARREF	T-Bar Tip Refill (3 tips)	€ 35,00
INLINEREF	In-Line Refill (3 tips)	€ 35,00
BLPGREF	Biological Periogauge (BLPG) Refill (3 tips)	€ 35,00
PAPREF	Papilla Refill (3 tips)	€ 35,00



ZVANI 8899 UN SAŅEM VAIRĀK PAR SAVU PASŪTĪJUMU!